

Self Seal Sterilisation Pouches



- Made in Britain by Granton Medical in a Class 10,000 Clean Room facility on unique, high speed, fully validated plant
- Heavy weight Medical Grade Paper
- Blue tinted transparent, multi-layer PET/PP film allows for easy viewing of the contents of the pouch
- Self Seal strip, quick and easy to use and eliminates the need for heat seal equipment and provides a validated seal
- Suitable for both steam (autoclave) and ethylene oxide (EO) sterilisation methods
- Each pouch printed with water based, non-toxic process indicators for both steam and ethylene oxide sterilisation, compliant with ISO 11140-1
- Each pouch printed with a lot number for full manufacturing traceability
- Off-Set webs for easy opening
- Tack seals to prevent curling and reduce external contamination
- Comprehensive size range, custom sizes available by arrangement
- Comply with ISO 11607 and EN 868-5

Standard sizes available from stock:-

CODE	NOMINAL SIZE	BOX QTY	CASE QUANTITY
11110	2¼" x 4" / 60mm x 100mm	Box 200	10 boxes of 200 per case
11120	3½" x 5½" / 90mm x 140mm	Box 200	10 boxes of 200 per case
11128	3½" x 8" / 90mm x 200mm	Box 200	10 boxes of 200 per case
11130	3½" x 9" / 90mm x 230mm	Box 200	10 boxes of 200 per case
11140	5¼" x 10" / 135mm x 255mm	Box 200	10 boxes of 200 per case
11144	5¼" x 14" / 135mm x 360mm	Box 200	10 boxes of 200 per case
11150	7½" x 12" / 190mm x 305mm	Box 200	5 boxes of 200 per case
11153	7½" x 13" / 190mm x 330mm	Box 200	5 boxes of 200 per case
11185	12" x 15" / 305mm x 381mm	Box 200	5 boxes of 200 per case
11190	5" x 15" / 130mm x 385mm	Pack 200	5 boxes of 200 per case
11193	8" x 15" / 200mm x 385mm	Pack 200	5 boxes of 200 per case
11198	10" x 17" / 250mm x 440mm	Pack 200	5 boxes of 200 per case
11199	12" x 17" / 300mm x 440mm	Pack 200	5 boxes of 200 per case



LOADING INSTRUMENT INTO POUCH



PEELING OFF STRIP FROM ADHESIVE BAND



FOLDING PAPER FLAP ALONG PERFORATION



PRESSING DOWN FIRMLY TO CREATE SEAL



PEELING OPEN FROM OPPOSITE END AFTER PROCESSING



Granton®
Medical